

Zoll AED Plus

A cardiac arrest or cardiac fibrillation is always sudden and, at that moment it is necessary to act immediately. Of course, one of the first things one should do is contact professional help. Zoll AED Plus does not only alert you when a transient is needed, also tells you if the heart massage is deep enough and runs at the right motion.



Features

For the victim however, it is of the highest importance that directly first aid is granted. This first aid exists in the base from CPR and rescue breathing; an AED can increase the survival rate of the victim.

Dräger has consciously chosen for the Zoll's AED Plus.

This AED distinguishes itself from other brands

- ignores vibrations and movements of the vessel, making a good analysis
- based on IP55
- makes use of special CPR-D pad; a pressure point mounted to the electrodes, with which the pressure and the pace of the heart massage is controlled
- a good step-by-step-resuscitation guidance with understandable icons and clear voice
- a lifetime of batteries and electrodes of 5 years without maintenance
- the Zoll AED Plus is also the only AED using normal batteries, these are easy to replace and cheap to run the Zoll AED has a standard LCD screen 7 x 3 cm with, optionally, text or ECG-view.

TECHNICAL SPECIFICATIONS

Particle and water ingress	IP-55
Battery life	5 years
Display	LCD display, 70 x 30 mm
Dimensions	133 x 241 x 292 mm
Weight	3.1 kg
Operating temperature	0°C to +50°C
Storage temperature	-30°C to +60°C
Cable length	1.1 m
Recharge time	Less than 10 seconds when using new batteries
Energy level shock	120J-150J-200J
Warranty	7 years
Approvals	UL2601, AAMI DF-39, IEC 601-2-4, EN 60601-1